

## Standard LensCheck Lens Adapters

- LensCheck Tower Mount (Carrier and Lens Mount Tower) + C-Mount Plate (1 inch x 32 TPI) Included with LensCheck System

### Additional Options Are:

- Blank Adapter Plate
- S-Mount - M12 x 0.5mm
- Self-Centering (Small Caliper) Lens Mount
- Canon EF-Mount
- D-mount 0.625"-32
- Leica M-Mount
- M13 Mount M13 X 1 mm
- M16 Mount M16 X 0.5mm
- M39 Mount M39 X 0.977mm
- M42 Mount M42 X 1mm
- Mitutoyo Mount M 26 X 0.706mm
- Nikon F-mount
- Nikon Mount M25 X 0.75mm
- Olympus OM-Mount
- Pentax K-Mount
- DIN Standard Mount 0.800" x 36 TPI (Formerly RMS)
- S Mount Lens Holder M12 X 0.5
- Sigma SA-Mount
- Thorlabs SM05 0.535"-40 TPI
- Thorlabs SM1 1.035"-40 TPI
- Thorlabs SM2 2.035"-40 TPI
- Thorlabs SM30 M30.5x0.5 TPI
- T-mount M42 X 0.75mm
- Zeiss Mount M27x0.75mm

*\*\*Is your lens mount not listed? Custom lens adapters are available on demand.*

## Additional LensCheck Mounts

### LensCheck Small Lens Mount

- Small opti-claw mount assembly with tip/tilt adjustment
- Compatible with existing LensCheck rail
- Appropriate for optics with diameters of 1.5 - 16 mm

### LensCheck Medium Lens Mount

- Small opti-claw mount assembly with tip/tilt adjustment
- Compatible with existing LensCheck rail
- Appropriate for optics with diameters up to 50 mm



The Optical Engineering Experts®

**Anywhere Light Goes®**

Optikos Corporation  
107 Audubon Rd., Bldg. 3,  
Wakefield, MA 01880 USA

**phone:** +1.617.354.7557  
**email:** sales@optikos.com  
**web:** optikos.com

*Copyright © 2023 Optikos Corporation; all rights reserved. Optikos®, OpTest®, Meridian®, VideoMTF®, The Optical Engineering Experts® and Anywhere Light Goes® are registered trademarks of Optikos Corporation.*